



SenTower Park

Name rider: _____

Nr boxes: _____

Instructions <400 horses:

1. Place your horses in the allocated barn. (with or without First bedding) Do not unload any tack or equipment
2. Put the Temperature monitor sheet (filled in) on each barn door and have your horses' FEI Self Certifying document (400CHF fine/horse if not filled in), passports and Vet declaration (free of Rhino symptoms) ready. (NO declaration = NO Stable entrance)
3. Bring this sheet back to the stable manager. The vet now knows he can come to your stable(s).
- 4) The vet scans the horse's microchip to confirm identity (using Bluetooth microchip reader and FEI HorseApp)
- 5) The vet checks for other signs of infectious diseases and vaccination against equine influenza, as already indicated by the Veterinary Regulations.
- 6) The groom checks the temperature of each horse using their own thermometer.
- 5) After checking, the stable / Temperature monitor sheet will receive a sticker:

Green (>36.5° / <38.5°) = Unloading tack and equipment

Orange (>38.5° / <39.0°)= On hold + recheck temperature after 1hr

Red = (>39.0°) Remove the horse and all in-contact horses to the quarantine stable